

# MICROPETTE PLUS SINGLE CHANNEL

## MECHANICAL AUTOCLAVABLE PIPETTE - MICROPETTE PLUS - SINGLE-CHANNEL ADJUSTABLE VOLUME

Codice	Descrizione	Volume
7030301001	 <p>MECHANICAL AUTOCLAVABLE PIPETTE - MICROPETTE PLUS SINGLE-CHANNEL ADJUSTABLE VOLUME</p>	0.1-2.5µl
7030301004		0.5-10µl
7030301005		2-20µl
7030301006		5-50µl
7030301008		10-100µl
7030301009		20-200µl
7030301011		50-200µl
7030301014		100-1000µl
7030301011		200-1000µl
7030301014		1000-5000µl

## MECHANICAL AUTOCLAVABLE PIPETTE - MICROPETTE PLUS - SINGLE-CHANNEL FIXED VOLUME

Codice	Descrizione	Volume
7030302018	 <p>MECHANICAL AUTOCLAVABLE PIPETTE - MICROPETTE PLUS SINGLE-CHANNEL FIXED VOLUME</p>	5µl
7030302019		10µl
7030302020		20µl
7030302021		25µl
7030302022		50µl
7030302023		100µl
7030302024		200µl
7030302025		250µl
7030302026		500µl
7030302027		1000µl
7030302028		2000µl
7030302029		5000µl

CON L'ACQUISTO  
DI ALMENO 4 PIPETTE,  
**IN OMAGGIO**  
**UNA MINI**  
**CENTRIFUGA**  
**DA BANCO**

**COD. 9031001012**



**D1008 MINI  
CENTRIFUGE  
WITH BLUE LID,  
INCLUDING TWO  
ROTORS  
EURO PLUG**

# MICROPETTE PLUS MULTI-CHANNEL



## MICROPETTE PLUS - MULTI-CHANNEL - EIGHT-CHANNEL ADJUSTABLE VOLUME

Codice	Descrizione	Volume
7030303004	MICROPETTE PLUS MULTI-CHANNEL EIGHT-CHANNEL ADJUSTABLE VOLUME	0.5-10µl
7030303006		5-50µl
7030303012		50-300µl



## MICROPETTE PLUS - MULTI-CHANNEL - TWELVE-CHANNEL ADJUSTABLE VOLUME

Codice	Descrizione	Volume
7030304004	MICROPETTE PLUS MULTI-CHANNE TWELVE-CHANNEL ADJUSTABLE VOLUME	0.5-10µl
7030304006		5-50µl
7030304012		50-300µl



Shipping  
Time  
**24h**



Shipping  
Fees  
**€ 0**



Shipping Fees for temperature  
controlled products  
**€ 0**



Minimun  
Order  
**NO**

# PIPETTATORE LEVO PLUS



## FEATURES

- Compatible with most plastic and glass pipettes from 0.1 to 100mL.
- Eight speeds of selection for aspiration and dispensing.
- LCD display shows low battery warning and speed setting.
- Powerful rechargeable Li-ion battery enables 8h runtime.
- Enables single handed operation with minimum effort 0.45µm replaceable hydrophobic filter.

## PIPETTE FILLER-LEVO PLUS

Codice	Descrizione
7033100100	LEVO PLUS WITH AC ADAPTER, EURO PLUG

## ACCESSORIES

Codice	Descrizione
17000205	LEVO PLUS LITHIUM BATTERY
17000103	FILTER 0.45UM, 1PC/PK

## ACCESSORIES

<b>ASPIRATION SPEED</b>	25ML<5S (8 SHIFT)
<b>DISPENSING SPEED</b>	MOTOR (8 SHIFT)/GRAVITY
<b>BATTERY</b>	LITHIUM-ION
<b>BATTERY SERVICE LIFE</b>	MORE THAN 8 HOURS OF INTERMITTENCE USE
<b>CHARGING TIME</b>	2-3 HOURS
<b>PIPETTE TYPES</b>	GLASS OR PLASTIC PIPETTE (0.1-100ML), PASTEUR PIPETTES
<b>FILTER</b>	0.45MM HYDROPHOBIC
<b>WEIGHT</b>	200G



Shipping  
Time  
**24h**



Shipping  
Fees  
**€ 0**



Shipping Fees for temperature  
controlled products  
**€ 0**



Minimun  
Order  
**NO**